

ILK antibody

Catalog Number: 111753

Product name

ILK antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ILK Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, MCF-7 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of ILK recombinant protein (Accession Number: NM_001014794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111753(ILK Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:111753(ILK Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)