

USP22 Polyclonal Antibody



Catalog Number: 164753

Product name

USP22 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to USP22

Preparation

Antigen: Recombinant protein of human USP22

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

Validations



Western blot - USP22 Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using USP22 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.