

# RNF213 antibody



Catalog Number: 114753

## Product name

RNF213 antibody

## Specificity

Human; other species not tested.

## Antibody description

RNF213 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human spleen tissue.

## Preparation

This antibody was obtained by immunization of RNF213 recombinant protein (Accession Number: BC032220). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human lymphoma slide using Catalog No:114753(RNF213 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human lymphoma slide using Catalog No:114753(RNF213 Antibody) at dilution of 1:50