

RNF213 antibody



Catalog Number: 114753

Product name

RNF213 antibody

Specificity

Human; other species not tested.

Antibody description

RNF213 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of RNF213 recombinant protein (Accession Number: BC032220). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human lymphoma slide using Catalog No:114753(RNF213 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human lymphoma slide using Catalog No:114753(RNF213 Antibody) at dilution of 1:50