

ACVR1B antibody



Catalog Number: 107753

Product name

ACVR1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACVR1B Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, DMSO treated HepG2 cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of ACVR1B recombinant protein (Accession Number: NM_004302). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107753(ACVR1B Antibody) at dilution of 1:300