

NDUFA7 antibody



Catalog Number: 113063

Product name

NDUFA7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDUFA7 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 13 kDa

Preparation

This antibody was obtained by immunization of NDUFA7 recombinant protein (Accession Number: NM_005001). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113063(NDUFA7 Antibody) at dilution of 1:300