

Anti-TrkC/NTRK3 antibody



Catalog Number: 104063

Product name

Anti-TrkC/NTRK3 antibody

Immunogen

[Human TrkC/NTRK3 \(His Tag\) recombinant protein](#)

Specificity

Human TrkC / NTRK3

Antibody description

Mouse monoclonal to TrkC/NTRK3

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human TrkC / NTRK3 (rh TrkC / NTRK3; NP_001007157.1; Met1-Asp428) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



TrkC / NTRK3 Antibody (PE), Mouse MAb, Flow cytometric

Profile of anti-TrkC reactivity on A549 cells analyzed by flow cytometry.