

Anti-14-3-3 alpha/beta Rabbit antibody



Catalog Number: 168063

Product name

Anti-14-3-3 alpha/beta Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

14-3-3 alpha/beta Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human 14-3-3 alpha/beta

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of 14-3-3 alpha/beta in HeLa, A549, MCF-7, U251, C2C12 cell lysates using 14-3-3 alpha/beta Rabbit mAb(1:1000 diluted). Predicted band size: 28kDa. Observed band size: 28kDa.