

# ANKMY2 antibody



Catalog Number: 108063

## Product name

ANKMY2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ANKMY2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of ANKMY2 recombinant protein (Accession Number: NM\_020319). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108063(ANKMY2 antibody) at dilution of 1:500