

Anti-P44/42 MAPK(ERK1/2) Mouse antibody



Catalog Number: 169063

Product name

Anti-P44/42 MAPK(ERK1/2) Mouse antibody

Dilutions

1:1000-2000

Specificity


1:50-100

Human, Rat, Mouse

Validations

Antibody description

P44/42 MAPK(ERK1/2) Mouse monoclonal antibody



Preparation

Antigen: Synthetic Peptide of P44/42 MAPK (ERK1/2)

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using P44/42 MAPK (ERK1/2) Mouse mAb diluted at 1:200.

Storage

-20°C/1 year

Clonality


Monoclonal

Ig Type

Mouse IgG

Applications

WB, IHC-P



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with P44/42 MAPK(ERK1/2) Mouse mAb diluted at 1:2,000.