

# CD1A antibody



Catalog Number: 109063

## Product name

CD1A antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CD1A Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, Jurkat cells, mouse liver tissue, mouse skin tissue, mouse thymus tissue, Raji cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of CD1A recombinant protein (Accession Number: BC031645). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109063(CD1A antibody) at dilution of 1:800