

HRH1 antibody

Catalog Number: 111463

Product name

HRH1 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HRH1 Rabbit Polyclonal antibody. Positive WB detected in human colon tissue, Jurkat cells. Positive IHC detected in human endometrial cancer tissue, human colon cancer tissue, human trachea tissue. Observed molecular weight by Western-blot: 56 kDa

Validations




Preparation

This antibody was obtained by immunization of HRH1 recombinant protein (Accession Number: NM_001098212). Purification method: Antigen affinity purified.

human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111463(HRH1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:111463(HRH1 antibody) at dilution of 1:100 (under 40x lens)



Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:111463(HRH1 antibody) at dilution of 1:100

HRH1 antibody



Catalog Number: 111463

(under 10x lens)