

Mouse CRABP2 (His Tag) recombinant protein



Catalog Number: 500863

General Information

Gene Name Synonym

Cellular retinoic acid-binding protein II

Protein Construction

A DNA sequence encoding the mouse CRABP2 (P22935) (Pro 2-Glu 138) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

> 96 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant mouse CRABP2 consisting of 148 amino acids and has a calculated molecular mass

of 17.1 kDa. rm CRABP2 migrates as an approximately 18 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CRABP2 / CRABPII Protein (His Tag) SDS-PAGE