

# ABO Polyclonal Antibody



Catalog Number: 162863

## Product name

ABO Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ABO

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 54-354 of human ABO (NP\_065202.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Western blot - ABO Polyclonal Antibody

Western blot analysis of extracts of K-562 cells, using ABO antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - ABO Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using ABO antibody . Blue: DAPI for nuclear staining.