

PHPT1 antibody



Catalog Number: 113863

Product name

PHPT1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PHPT1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, human kidney tissue. Observed molecular weight by Western-blot: 14kd

Preparation

This antibody was obtained by immunization of PHPT1 recombinant protein (Accession Number: NM_014172). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113863(PHPT1 antibody) at dilution of 1:500