

TBC1D19 antibody



Catalog Number: 115863

Product name

TBC1D19 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBC1D19 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human pancreas tissue, human testis tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of TBC1D19 recombinant protein (Accession Number: NM_018317). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IF, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115863(TBC1D19 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:115863(TBC1D19 Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded

TBC1D19 antibody



Catalog Number: 115863

human pancreas tissue slide using Catalog No:115863(TBC1D19 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:115863(TBC1D19 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)