

# Anti-NDUFS6 Rabbit antibody



Catalog Number: 166863

## Product name

Anti-NDUFS6 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

NDUFS6 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human NDUFS6

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/200

## Validations



Western blot detection of NDUFS6 in COS7,Raji,3T3,Jurkat cell lysates using NDUFS6 Rabbit pAb(1:1000 diluted).Predicted band size:14KDa.Observed band size:14KDa.