

Anti-CD22 antibody



Catalog Number: 142163

Product name

Anti-CD22 antibody

Immunogen

[Mouse CD22 \(His Tag\) recombinant protein](#)

Specificity

Mouse CD22

Antibody description

Rabbit Monoclonal to Mouse CD22

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD22 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 0.5-2 µg/Test

Validations



CD22 Antibody, Rabbit MAb, Flow Cytometry