

ZNF34 antibody

Catalog Number: 117163

Product name

ZNF34 antibody

Specificity

human; other species not tested.

Antibody description

ZNF34 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, Jurkat cells. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 62-68kd

Preparation

This antibody was obtained by immunization of ZNF34 recombinant protein (Accession Number: XM_011517318). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117163(ZNF34 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:117163(ZNF34 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:117163(ZNF34 antibody) at dilution of 1:50 (under 40x lens)