

# ZNF34 antibody

Catalog Number: 117163

## Product name

ZNF34 antibody

## Specificity

human; other species not tested.

## Antibody description

ZNF34 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, Jurkat cells. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 62-68kd

## Preparation

This antibody was obtained by immunization of ZNF34 recombinant protein (Accession Number: XM\_011517318). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117163(ZNF34 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:117163(ZNF34 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:117163(ZNF34 antibody) at dilution of 1:50 (under 40x lens)