

# Anti-Cytokeratin 19 (10G10) Mouse antibody



Catalog Number: 168163

## Product name

WB, IHC-P

Anti-Cytokeratin 19 (10G10) Mouse antibody

## Dilutions

## Specificity

IHC: 1/200

Human, Mouse, Rat

## Antibody description

Cytokeratin 19 (10G10) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic Peptide of CK19

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Mouse IgG1

## Applications

## Validations



Western blot analysis of extracts from PC-3, SKOV-3 and MCF7 cells using Keratin 19 Mouse mAb at 1:1000 dilution. Predicted band size: 44kDa. Observed band size: 44kDa.



IHC staining of human colon cancer tissue with CK19 mouse mAb(11F5) diluted at 1:200.