

TCEAL8 Polyclonal Antibody



Catalog Number: 165563

Product name

TCEAL8 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TCEAL8

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human TCEAL8 (NP_699164.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100

Validations



Immunofluorescence - TCEAL8 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TCEAL8 antibody . Blue: DAPI for nuclear staining.