

CRISP2 antibody



Catalog Number: 109563

Product name

CRISP2 antibody

Specificity

human, mouse; other species not tested.

Antibody description

CRISP2 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Positive IHC detected in human testis. Observed molecular weight by Western-blot: 27kd

Preparation

This antibody was obtained by immunization of CRISP2 recombinant protein (Accession Number: NM_003296). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109563(CRISP2 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:109563(CRISP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:109563(CRISP2 antibody) at dilution of 1:50 (under 40x lens)