

# KCNG4 antibody

Catalog Number: 111963

## Product name

KCNG4 antibody

## Specificity

Human; other species not tested.

## Antibody description

KCNG4 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of KCNG4 recombinant protein (Accession Number: NM\_133490). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111963(KCNG4 Antibody) at dilution of 1:700