

Anti-CD81 Rabbit antibody



Catalog Number: 166963

Product name

Anti-CD81 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CD81 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human CD81

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC


Dilutions

WB: 1/1000-1/2000

Validations



Immunohistochemical analysis of paraffin-embedded human tonsil, using CD81 Antibody.



Western blot analysis of CD81 expression in Ramos cell lysate.