

PTCH2 antibody

Catalog Number: 114263

Product name

PTCH2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTCH2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human brain tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 131-140 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003738). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114263(PTCH2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:114263(PTCH2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114263(PTCH2 antibody) at dilution of 1:50 (under 40x lens)