

SHANK1 antibody



Catalog Number: 115263

Product name

SHANK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SHANK1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 250-300 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_016148). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115263(SHANK1 antibody) at dilution of 1:300