

Anti-POLR3H antibody



Catalog Number: 176263

Product name

Anti-POLR3H antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to POLR3H

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rpc8. aa range:151-200

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of K562 cells using POLR3H Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).