

# Anti-Ran Rabbit antibody



Catalog Number: 167263

## Product name

Anti-Ran Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Ran Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Ran

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

FC: 1/20-1/100

IP: 1/20-1/50

## Validations



Western blot detection of Ran in Jurkat, HeLa, COS7, Raji cell lysates using Ran Rabbit pAb(1:1000 diluted). Predicted band size: 24KDa. Observed band size: 24KDa.