

ZPBP antibody

Catalog Number: 117263

Product name

ZPBP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZPBP Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of ZPBP recombinant protein (Accession Number: NM_007009). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-ZPBP (IP: Catalog No: 117263, 3ug; Detection: Catalog No: 117263 1:300) with mouse testis tissue lysate 2800ug.