

PDE4D antibody



Catalog Number: 113663

Product name

PDE4D antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDE4D Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse pancreas tissue. Positive IHC detected in mouse colon tissue. Observed molecular weight by Western-blot: 91 kDa

Preparation

This antibody was obtained by immunization of PDE4D recombinant protein (Accession Number: NM_001197223). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:50-1:500

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113663(PDE4D antibody) at dilution of 1:400