

Anti-Phospho-CDK2 (Thr160) Rabbit antibody



Catalog Number: 166663

Product name

Anti-Phospho-CDK2 (Thr160) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Phospho-CDK2 (Thr160) Rabbit polyclonal antibody

Dilutions

WB: 1/500

Preparation

Antigen: CVPVRTYp-THEVVTL

Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Validations



Western blot analysis of extracts from HeLa cells, untreated or UV-treated, using Phospho-CDK2 (Thr160) Rabbit pAb (1:500 diluted). Predicted band size: 33kDa. Observed band size: 33kDa..