

SLC27A3 antibody

Catalog Number: 115363

Product name

SLC27A3 antibody

Immunogen

[Human SLC27A3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC27A3 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in human testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 79 kDa

Preparation

This antibody was obtained by immunization of SLC27A3 recombinant protein (Accession Number: BC029792). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115363(SLC27A3 antibody) at dilution of 1:500



IP Result of anti-SLC27A3 (IP:Catalog No:115363, 4ug; Detection:Catalog No:115363 1:500) with mouse testis tissue lysate 4000ug.