

Human UBE2H (GST tag) recombinant protein



Catalog Number: 519363

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 H, UbcH2, Ubiquitin Carrier Protein H, Ubiquitin-Conjugating Enzyme E2-20K, Ubiquitin-Protein Ligase H, UBE2H

Protein Construction

A DNA sequence encoding the Human UBE2H Met1-Leu183 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2H (GST tag) recombinant protein