

# ZNF707 Polyclonal Antibody



Catalog Number: 161763

## Product name

ZNF707 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ZNF707

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-180 of human ZNF707 (NP\_776192.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Immunohistochemistry - ZNF707 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using ZNF707 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ZNF707 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using ZNF707 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZNF707 Polyclonal Antibody

# ZNF707 Polyclonal Antibody



Catalog Number: 161763

---

Immunohistochemistry of paraffin-embedded

human gastric cancer using ZNF707 antibody at dilution of 1:100 (40x lens).