

SYN1-Specific antibody

Catalog Number: 115763

Product name

SYN1-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SYN1-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in rat brain tissue. Observed molecular weight by Western-blot: 70-80 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006950). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IP: 1:1000-1:10000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115763(SYN1-Specific Antibody) at dilution of



1:5000

IP Result of anti-SYN1-Specific (IP:Catalog No:115763, 4ug; Detection:Catalog No:115763 1:3000) with rat brain tissue lysate 4000ug.