

# VMAT2 antibody



Catalog Number: 116763

## Product name

VMAT2 antibody

## Immunogen

[Human VMAT2 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

VMAT2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, SGC-7901 cells. Observed molecular weight by Western-blot: 65-75kd,55-60kd

## Preparation

This antibody was obtained by immunization of VMAT2 recombinant protein (Accession Number: NM\_003054). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116763(VMAT2 antibody) at dilution of 1:300