

VMAT2 antibody



Catalog Number: 116763

Product name

VMAT2 antibody

Clonality

Polyclonal

Immunogen

[Human VMAT2 Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

VMAT2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, SGC-7901 cells. Observed molecular weight by Western-blot: 65-75kd, 55-60kd


Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of VMAT2 recombinant protein (Accession Number: NM_003054). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116763(VMAT2 antibody) at dilution of 1:300