

Anti-N Cadherin Rabbit antibody



Catalog Number: 167763

Product name

Anti-N Cadherin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

N Cadherin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human N Cadherin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC: 1/20

Validations



Western blot detection of N Cadherin in HeLa,A549,U2OS,U87-MG cell lysates using N Cadherin Rabbit mAb(1:1000 diluted).Predicted band size:100kDa.Observed band size:140kDa.