

CACNB2 antibody



Catalog Number: 108763

Product name

CACNB2 antibody

Immunogen

[Human CACNB2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

CACNB2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of CACNB2 recombinant protein (Accession Number: NM_201596). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108763(CACNB2 Antibody) at dilution of 1:300