

# NDUFB3 antibody



Catalog Number: 113073

## Product name

NDUFB3 antibody

## Immunogen

[Human NDUFB3 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NDUFB3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 11 kDa

## Preparation

This antibody was obtained by immunization of NDUFB3 recombinant protein (Accession Number: NM\_001257102). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113073(NDUFB3 antibody) at dilution of 1:500