

NDUFB3 antibody



Catalog Number: 113073

Product name

NDUFB3 antibody

Clonality

Polyclonal

Immunogen

[Human NDUFB3 Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

NDUFB3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 11 kDa

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of NDUFB3 recombinant protein (Accession Number: NM_001257102). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113073(NDUFB3 antibody) at dilution of 1:500

