

CKS1B Polyclonal Antibody



Catalog Number: 165073

Product name

CKS1B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CKS1B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human CKS1B (NP_001817.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - CKS1B Polyclonal Antibody

Western Blot analysis of Jurkat cell nucleus extracts using Cks1 Polyclonal Antibody.