

# Anti-C1orf35 antibody

Catalog Number: 176073

---

## Product name

Anti-C1orf35 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to C1orf35

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mmtag2. aa range:10-59

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.;; Western blot analysis of 293T lysis using MMTAG2 antibody. Antibody was diluted at 1:1000