

HMBS antibody



Catalog Number: 111473

Product name

HMBS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HMBS Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A375 cells, HEK-293 cells, Raji cells. Positive IF detected in MCF-7 cells, Hela cells, HepG2 cells. Observed molecular weight by Western-blot: 39-42kd

Preparation

This antibody was obtained by immunization of HMBS recombinant protein (Accession Number: NM_000190). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111473(HMBS antibody) at dilution of 1:300



Immunofluorescent analysis of MCF-7 cells, using HMBS antibody Catalog No:111473 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).