

# CLDN17 Polyclonal Antibody



Catalog Number: 162473

## Product name

CLDN17 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to CLDN17

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human CLDN17 (NP\_036263.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Validations



Western blot - CLDN17 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CLDN17 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.