

# APOBEC3F Polyclonal Antibody



Catalog Number: 163473

## Product name

APOBEC3F Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to APOBEC3F

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 189-373 of human APOBEC3F (NP\_660341.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - APOBEC3F Polyclonal Antibody

Western blot analysis of extracts of PC3 cells, using APOBEC3F antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.