

# PAAF1 antibody

Catalog Number: 113473

## Product name

PAAF1 antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

## Specificity


Human, Mouse; other species not tested.

IF: 1:20-1:200

## Antibody description

PAAF1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells, human brain tissue, Jurkat cells. Positive IHC detected in human testis tissue, human brain tissue, human kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

## Validations




## Preparation

This antibody was obtained by immunization of PAAF1 recombinant protein (Accession Number: NM\_001267805). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113473(PAAF1 antibody) at dilution of 1:2000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113473(PAAF1 Antibody) at dilution of 1:50 (under 10x lens)




# PAAF1 antibody

Catalog Number: 113473

---

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113473(PAAF1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HEK-293 cells using Catalog No:113473(PAAF1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)