

SNX18 antibody



Catalog Number: 115473

Product name

SNX18 antibody

Immunogen

[Human SNX18 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SNX18 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, human placenta tissue, L02 cells, MCF7 cells. Observed molecular weight by Western-blot: 67-70kd

Preparation

This antibody was obtained by immunization of SNX18 recombinant protein (Accession Number: NM_001102575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115473(SNX18 antibody) at dilution of 1:1000