

# Anti-CLEC4A antibody



Catalog Number: 176473

## Product name

Anti-CLEC4A antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to CLEC4A

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human clec4a. aa range:161-210

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:: Western Blot analysis of HEPG2 cells using CLECSF6 Polyclonal Antibody diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human colon cancer, antibody was diluted at 1:200